

## Protein Dimerization and Aggregation

Unexpected high molecular weights on western blots have observed for many proteins. They are may be caused by dimerization and aggregation.

### Preparation of Whole Cell Lysate

Lysis buffer:

1% Triton X-100

0.5% Nonidet P-40

150mM NaCl (58.44g/mol)

10mM Tris (121.14g/mol), pH 7.4

2mM EDTA (372.2g/mol)

1mM EGTA (380.4g/mol)

Protease inhibitors (83µg/ml aprotinin, 30 µg/ml leupeptin, 1mg/ml Pefabloc, 50µg/ml calpain inhibitor, 50µg/ml bestatin and 5µg/ml pepstatin.

It is also important to add 100mM iodoacetamide (Sigma I1149-5G, MWT 184.96g/mol stored at 4C) to the lysis buffer to block any non-specific sulfide link that might form between free cystines during the lysis and extraction of the proteins (i.e. avoid aggregate formation). Iodoacetamide will block thiol of Cys , thus, if aggregate is still appears after iodacetamide addition to the lysis buffer, this means that the cys-cys bond is SPECIFIC.

Procedure:

- (1) Prepare the lysis buffer with and without 100mM Iodoacetamide (to examine the non specific disulfide link).
- (2) Cells are washed twice with PBS
- (3) Cells are solubilized with the lysis buffer O/N at 4C using end-over-end rotor.
- (4) Spin the lysate at 15,000rpm for 15min at 4C (to remove insoluble stuff).
- (5) The supernatant is considered as a total cell lysate
- (6) Determine the protein concentration in the sample (use Lowery method).

Thus, to prepare 50ml Lysis buffer (don't add iodoacetamide)

500µl Triton X-100

250µl NP-40

0.438g NaCl

0.06g Tris (pH 7.4)

0.0372g EDTA

0.019g EGTA

Add the protease inhibitor cocktail (Sigma P-8340) as 1:100 (i.e. 10µl of the protease for each 1ml of the lysate). Add 0.0925g of Iodoacetamide for 5ml of the lysis buffer to get 0.1M .